

Catalog Number: 230-00523

Recombinant Alt a1 allergen protein**Source**


Species	Alternaria alternata
Accession Number	Q6Q128
Gene Symbol	ALT A1
Expressed Region	Gln2-Ser157
Synonyms	Alt a1

Preparation

Expression System	E.Coli
Tag	N-terminal 6x histidine tag
Purification	His-tag affinity purification by immobilized metal ion affinity chromatography (IMAC)
Purity	> 95%
Purity Determined By	SDS-PAGE under reducing conditions and visualized by Coomassie blue staining
Molecular Weight	16.9

Protein Specifications

Format	Liquid
Formulation	Supplied as a 0.2 um filtered solution in 2 M urea, 50 mM Tris-HCl (pH 8.0).
Concentration	Lot specific (see the label on the vial)
Preservative	None
Endotoxin Level	Not determined



Shipping

Ice packs

Storage/Stability

Upon arrival, the protein may be stored for 2 weeks at 4 °C. For long term storage, it is recommended to store at -20 °C or -80 °C in appropriate aliquots. Avoid repeated freeze-thaw cycles.

This product is furnished for LABORATORY RESEARCH USE ONLY.

Not for diagnostic or therapeutic use.